Wolfram Droh GmbH

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

MD

REV GA /2021-04

Ningbo Huakun Medical Equipment Co., Ltd.
No. 208 Binjiang Road
Cixi High-tech Industrial Development Zone
Cixi City
315301 Zhejiang, P.R. China

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GB Contents

1	General	Page	11
2	Safety and Responsibilities	Page	11
3	Modification of the product	Page	12
4	Before using the product	Page	12
5	Using the product	Page	13
6	Product Maintenance / Reprocessing	Page	13
7	Accessories	Page	13
8	Replacement parts and repairs	Page	14
9	Disposal	Page	14
10	Technical Specifications	Page	14
11	Explanation of symbols	Page	15

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1 General

- This product may only be used by trained professionals under medical instruction.
- The catheter mount is intended to be used for the transportation of medical gases to patient and to enable a better mobility of the patient. Usually the product is used as a connection between the intubator and the breathing system.
- Some versions have suction openings to allow suctioning from the patient using a bronchial catheter during intubation.
- It is mainly used in hospitals in the operating room for anaesthesia and machine ventilation.
- These instructions reflect the state of the art at the time of publication.
- This product complies with the provisions of the EU Directive 93/42/EEC (Medical Device Directive) and Medical Device Regulation (EU) 2017/745.
- The instructions for use comply with DIN EN 82079-1.
- This product meets the requirements of DIN EN 5367.

2 Safety and Responsibilities

IMPORTANT! READ CAREFULLY BEFORE USE

KEEP FOR FUTURE REFERENCE

- Please observe the following instructions and care instructions before use.
- If in doubt about the use and connection of the product, please contact customer service.
- If additional oxygen is administered care should be taken to ensure that smoking and other naked flames are forbidden in the vicinity of oxygen-carrying equipment and fittings.
- No flammable anaesthetic gases may be used due to possible static charges.
- Keep the produce free of oil and grease, otherwise there is a risk of explosion. This includes the exposure to cleaning alcohols, hand cream and adhesives – in contact with highly compressed oxygen explosive reactions may occur.
- The product must not come into contact with a pointed or a sharp object which could damage it. This could lead to breakages or similar damage and the product would start to leak.
- The product may not be used for ventilating babies.
- Make sure that the correct size is used for each individual patient.

- **Single use products may not be reused!**
They must be disposed of immediately after use in accordance with chapter 9.
The reuse of single use products represents a potential risk for the patient or the user. It can result in contamination and/or reduced functionality.
- Contamination and/or limited functionality of the product can lead to injury, illness or death of the patient.

3 Modification of the product

- Apart from the described maintenance instructions, the product may not be modified or repaired by the operator!
- The manufacturer's responsibility is limited to the dimensional and physiological suitability of the tubing.
- Improper use, modifications to the product or negligence will void any warranty. This also applies to unauthorised repairs. Normal signs of wear and tear on the device are also excluded.

4 Before using the product

- Do not expose the product to direct sunlight or other light sources. They must be stored and transported in a dry place and in their original packaging.
- Remove all delivered parts from the packaging. If necessary, remove the respective plastic covers and dispose of the packaging according to the current national regulations.
- Before each use verify that parts or components do not show any signs of hardening, damage, breakage or cracks. The product must also not be contaminated or soiled. If any of the conditions described above are present, the device must not be put into operation and must be disposed of in accordance with chapter 9.
- In case of doubt perform a comparison check with a brand-new product.
→ **Damaged products must not be used!**
They must be disposed of in accordance with the instructions in chapter 9.
- The product contains PP, PE and EVA.
- Please note that the product is supplied in a non-sterile condition.
- The hygienic preparation of the product lies within the clinical responsibility of the user! The plastic product can be disinfected with a suitable cleaning agent (e.g. alcohol) before use.
Sterilization is not permitted!

- To attach and detach to/from the system, the product must be held at the connector to prevent cracks from forming between the connector and the main tubing.
- The ventilation and anaesthesia hoses may only be connected to compatible connectors (e.g. according to DIN ISO 5356 Part 1).
- After assembly, the tubing must be checked for sufficient pull-off force and for kinks as well as any cracks or hardening.
- A leak test must be carried out before each use:
 - Connect the terminal end of connection hose the to the system.
 - Seal the other end.
 - Observe whether the pressure of 60 cmH₂O drops or air escapes within one minute.
 - If this is the case, the device must not be put into operation and must be disposed of in accordance with chapter 9.
- Furthermore, the product must be checked for foreign bodies and lumen occlusion.

5 Using the product

- Recommended maximum working temperature is 45°C.
- Recommended maximum working pressure is 35 kPa.
- After connection to a suitable terminal device, the product is ready for use.
- The product is suitable for use with oxygen as well as for air, anaesthetic gases and their gas mixtures.
- There must be no interaction of the gases with the product materials. Please refer to the product materials.

6 Product maintenance/Reprocessing

- Single use products may not be reused!
- The product must not be used on more than one patient.
- After termination of the application, the product must be disposed of according to chapter 9.

7 Accessories

- The user is responsible for the correct combination with other medical products.

8 Replacement parts and repairs

- Repairs must only be carried out by the manufacturer or by authorised persons.
- In case of functional or material problems, contact customer service immediately.
- If the cause of the error cannot be eliminated, the product must not be put into operation. After consulting the customer service, the product must be disposed of according to chapter 9 or returned to the sales distributor.

9 Disposal

- The product must be disposed of after use.
- Please note the expiry date if applicable.
The product must not be used after the expiry date!
- Safe disposal is carried out in accordance with current national regulations

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